

血细胞分析仪用校准物参考值表

型号：C3

批号：2503C

生产日期：2025.04.07

有效期至：2025.05.14

仪器	参数	单位	靶值	允许偏差
RT-7200	WBC	$\times 10^9/L$	10.7	± 0.3
	RBC	$\times 10^{12}/L$	4.31	± 0.1
	HGB	g/L	140	± 2
	MCV	fL	91.7	± 2
	PLT	$\times 10^9/L$	220	± 15
RT-7300	WBC	$\times 10^9/L$	11.2	± 0.3
	RBC	$\times 10^{12}/L$	4.28	± 0.1
	HGB	g/L	141	± 2
	MCV	fL	90.2	± 2
	PLT	$\times 10^9/L$	237	± 15
RT-7600 RT-7600S	WBC	$\times 10^9/L$	10.3	± 0.3
	RBC	$\times 10^{12}/L$	4.28	± 0.1
	HGB	g/L	137	± 2
	MCV	fL	89.7	± 2
	PLT	$\times 10^9/L$	241	± 15
RT-7900	WBC	$\times 10^9/L$	10.6	± 0.3
	RBC	$\times 10^{12}/L$	4.29	± 0.1
	HGB	g/L	136	± 2
	MCV	fL	86.6	± 2
	PLT	$\times 10^9/L$	228	± 15
RT-7600Vet	WBC	$\times 10^9/L$	10.4	± 0.3
	RBC	$\times 10^{12}/L$	4.04	± 0.1
	HGB	g/L	119	± 2
	MCV	fL	92.6	± 2
	PLT	$\times 10^9/L$	168	± 15