

血细胞分析仪用校准物参考值表

型号：C3

批号：2504C

生产日期：2025-5-6

有效期至：2025-6-19

仪器	参数	单位	靶值	允许偏差
RT-7200	WBC	$\times 10^9/L$	9.7	± 0.3
	RBC	$\times 10^{12}/L$	4.25	± 0.1
	HGB	g/L	134	± 2
	MCV	fL	86.9	± 2
	PLT	$\times 10^9/L$	212	± 15
RT-7300	WBC	$\times 10^9/L$	10.6	± 0.3
	RBC	$\times 10^{12}/L$	4.39	± 0.1
	HGB	g/L	140	± 2
	MCV	fL	87.4	± 2
	PLT	$\times 10^9/L$	228	± 15
RT-7600 RT-7600S	WBC	$\times 10^9/L$	9.5	± 0.3
	RBC	$\times 10^{12}/L$	4.24	± 0.1
	HGB	g/L	133	± 2
	MCV	fL	87.7	± 2
	PLT	$\times 10^9/L$	227	± 15
RT-7900	WBC	$\times 10^9/L$	9.8	± 0.3
	RBC	$\times 10^{12}/L$	4.25	± 0.1
	HGB	g/L	132	± 2
	MCV	fL	83.6	± 2
	PLT	$\times 10^9/L$	210	± 15
RT-7600Vet	WBC	$\times 10^9/L$	9.7	± 0.3
	RBC	$\times 10^{12}/L$	4.04	± 0.1
	HGB	g/L	116	± 2
	MCV	fL	90.6	± 2
	PLT	$\times 10^9/L$	161	± 15