

血细胞分析仪用质控物（光学法）参考值表

型号：RT-D5

批号：52503L-N-H

生产日期：2025.05.19

有效期至：2025.08.16

仪器	参数	单位	低值 (Lot.: 52503L)		中值 (Lot.: 52503N)		高值 (Lot.: 52503H)	
			参考值	偏差范围	参考值	偏差范围	参考值	偏差范围
Hemaray 50CRP、 Hemaray 51、 Hemaray 52、 Hemaray 53、 Hemaray 54、 Hemaray 55CRP、 Hemaray 56CRP、 Hemaray 57CRP、 Hemaray 58CRP、 Hemaray 59CRP	WBC	×10 ⁹ /L	3.7	±0.5	8.4	±1.0	19.9	±2.5
	NEU	%	47.5	±15	54.5	±15	61.3	±10
	LYM	%	35.4	±15	27.4	±15	17.6	±10
	MON	%	9.1	±5	8.4	±5	8.5	±5
	EOS	%	6.9	±5	8.7	±5	11.3	±5
	BAS	%	1.1	±1.1	1.0	±1.0	1.3	±1.3
	NEU	×10 ⁹ /L	1.76	±0.5	4.58	±1.5	12.21	±2.8
	LYM	×10 ⁹ /L	1.31	±0.5	2.30	±1.5	3.50	±2.8
	MON	×10 ⁹ /L	0.34	±0.2	0.71	±0.5	1.70	±1.0
	EOS	×10 ⁹ /L	0.26	±0.2	0.73	±0.5	2.25	±1
	BAS	×10 ⁹ /L	0.04	±0.04	0.09	±0.09	0.25	±0.25
	RBC	×10 ¹² /L	2.29	±0.18	4.74	±0.24	5.46	±0.30
	HGB	g/L	62	±4	135	±6	170	±8
	MCHC	g/L	348	±40	350	±40	359	±40
	MCH	pg	27.1	±3.0	28.5	±3.0	31.1	±3.0
	MCV	fL	77.7	±5	81.4	±5	86.6	±5
	RDW-CV	%	13.3	±3.0	13.8	±3.0	14.1	±3.0
	HCT	%	17.8	±2.5	38.6	±3.5	47.3	±4.5
	PLT	×10 ⁹ /L	50	±20	232	±40	503	±60
	MPV	fL	9.1	±3.0	9.3	±3.0	9.2	±3.0
Hemaray 83 Hemaray 86 Hemaray 89 Hemacounter SL	WBC	×10 ⁹ /L	3.8	±0.5	8.7	±1.0	21.2	±2.5
	NEU	%	49.1	±15	55.2	±15	61.4	±10
	LYM	%	37.5	±15	30.0	±15	20.5	±10
	MON	%	6.7	±5	6.5	±5	6.5	±5
	EOS	%	5.9	±5	7.5	±5	10.8	±5
	BAS	%	0.8	±0.8	0.8	±0.8	0.8	±0.8
	NEU	×10 ⁹ /L	1.90	±0.5	4.89	±1.5	13.17	±2.8
	LYM	×10 ⁹ /L	1.43	±0.5	2.62	±1.5	4.34	±2.8
	MON	×10 ⁹ /L	0.25	±0.2	0.57	±0.5	1.37	±1.0
	EOS	×10 ⁹ /L	0.22	±0.2	0.65	±0.5	2.28	±1
	BAS	×10 ⁹ /L	0.02	±0.02	0.01	±0.01	0.02	±0.02
	RBC	×10 ¹² /L	2.19	±0.18	4.46	±0.24	5.13	±0.30
	HGB	g/L	61	±4	131	±6	162	±8
	MCHC	g/L	335	±40	339	±40	339	±40
	MCH	pg	27.9	±3.0	29.4	±3.0	31.6	±3.0
	MCV	fL	83.0	±5	86.7	±5	93.2	±5
	RDW-CV	%	16.7	±3.0	15.5	±3.0	15.4	±3.0
	HCT	%	18.2	±2.5	38.7	±3.5	47.8	±4.5
	PLT	×10 ⁹ /L	48	±20	242	±40	567	±60
	MPV	fL	8.2	±3.0	8.1	±3.0	8.0	±3.0