



生化试剂出厂检验报告

Chemistry Reagent Inspection Report

产品名称 Product Name		腺苷脱氨酶测定试剂盒（酶法） Adenosine Deaminase Assay Reagent Kit (Enzymatic)		
规格 Specification		R1: 45ml×1, R2: 15ml×1, 校准: 1ml×1	产品型号 Product Model	/
批号 Batch No.		20250303	有效期 Expiration Date	2026/3/2
生产日期 Production Date		2025. 03. 03	检验日期 Inspection Date	2025. 03. 13
试剂批号 Reagent batch number		20250303	试剂有效期 Expiration date of Reagent	2026/3/2
校准/质控批号 Calibrator/Control batch number		20250303	校准有效期 Expiration date of Calibrator	2026/9/2
检验项目 Inspection Items		检验标准 Inspection Standard		检验结果 Inspection Results
判定 Judgment				
1	外观 Appearance	R1无色、R2无色或淡黄色的澄清液体，无沉淀及絮状物 R1 is colorless, clear liquid and R2 is colorless or yellowish, clear liquid, without sediments or floccules.		符合 Meet with the requirement
2	净含量 Net volume	液体试剂的净含量应不少于标示值 The net content of liquid reagents should not be less than the indicated value.		符合 Meet with the requirement
3	准确度 Accuracy	相对偏差不超过±10% The relative deviation should not exceed ±10%.		1.62% Pass
4	重复性 Repeatability	变异系数CV≤5% Coefficient of Variation (CV) ≤5%		1.17% Pass
5	线性范围 Linearity range	在0U/L~185U/L (37℃) 范围内，相关系数(r)应≥0.990 Within the range of 0U/L ~ 185U/L (37℃), the correlation coefficient (r) should be ≥0.990.		0.9996 Pass
结论 Conclusion		合格 Pass		
检验员 Inspector		莫宗儒	复核人 Reviewer	高曼姿
日期 Date		2025. 03. 13	日期 Date	2025. 03. 13

注：若无内容填写应用“/”。

Note: If there is no information, please fill with "/".

表格编号/Form Number: RT-RFM-Q8503

REV: A. 2/ 保存期限: 10年/Storage Life: 10 years

