

血细胞分析仪用质控物（光学法）参考值表

型号：RT-D5

批号：52504L-N-H

生产日期：2025.07.07

有效期至：2025.10.04

仪器	参数	单位	低值 (Lot.: 52504L)		中值 (Lot.: 52504N)		高值 (Lot.: 52504H)	
			参考值	偏差范围	参考值	偏差范围	参考值	偏差范围
Hemaray 50CRP、 Hemaray 51、 Hemaray 52、 Hemaray 53、 Hemaray 54、 Hemaray 55CRP、 Hemaray 56CRP、 Hemaray 57CRP、 Hemaray 58CRP、 Hemaray 59CRP	WBC	×10 ⁹ /L	4.0	±0.5	9.3	±1.0	20.5	±2.5
	NEU	%	47.0	±15	52.6	±15	59.5	±10
	LYM	%	33.8	±15	28.2	±15	18.2	±10
	MON	%	10.9	±5	9.2	±5	10.1	±5
	EOS	%	6.0	±5	8.5	±5	10.7	±5
	BAS	%	2.3	±2.3	1.5	±1.5	1.5	±1.5
	NEU	×10 ⁹ /L	1.88	±0.5	4.89	±1.5	12.20	±2.8
	LYM	×10 ⁹ /L	1.37	±0.5	2.61	±1.5	3.73	±2.8
	MON	×10 ⁹ /L	0.44	±0.2	0.85	±0.5	2.06	±1.0
	EOS	×10 ⁹ /L	0.24	±0.2	0.79	±0.5	2.20	±1
	BAS	×10 ⁹ /L	0.09	±0.09	0.14	±0.14	0.31	±0.31
	RBC	×10 ¹² /L	2.39	±0.18	4.69	±0.24	5.40	±0.30
	HGB	g/L	63	±4	136	±6	168	±8
	MCHC	g/L	344	±40	343	±40	346	±40
	MCH	pg	26.4	±3.0	29.0	±3.0	31.1	±3.0
	MCV	fL	76.6	±5	84.6	±5	89.8	±5
	RDW-CV	%	14.1	±3.0	13.1	±3.0	13.4	±3.0
	HCT	%	18.3	±2.5	39.7	±3.5	48.5	±4.5
	PLT	×10 ⁹ /L	49	±20	240	±40	515	±60
	MPV	fL	9.3	±3.0	9.3	±3.0	9.4	±3.0
Hemaray 83 Hemaray 86 Hemaray 89 Hemacounter SL	WBC	×10 ⁹ /L	3.8	±0.5	9.2	±1.0	21.9	±2.5
	NEU	%	50.3	±15	54.2	±15	57.5	±10
	LYM	%	38.6	±15	32.9	±15	22.3	±10
	MON	%	5.1	±5	4.0	±4	9.7	±5
	EOS	%	5.6	±5	8.5	±5	10.1	±5
	BAS	%	0.4	±0.4	0.4	±0.4	0.4	±0.4
	NEU	×10 ⁹ /L	1.91	±0.5	4.99	±1.5	12.59	±2.8
	LYM	×10 ⁹ /L	1.47	±0.5	3.02	±1.5	4.87	±2.8
	MON	×10 ⁹ /L	0.19	±0.19	0.36	±0.36	2.12	±1.0
	EOS	×10 ⁹ /L	0.21	±0.2	0.78	±0.5	2.21	±1
	BAS	×10 ⁹ /L	0.02	±0.02	0.04	±0.04	0.09	±0.09
	RBC	×10 ¹² /L	2.23	±0.18	4.33	±0.24	5.03	±0.30
	HGB	g/L	59	±4	130	±6	160	±8
	MCHC	g/L	335	±40	349	±40	349	±40
	MCH	pg	26.5	±3.0	30.0	±3.0	31.8	±3.0
	MCV	fL	78.9	±5	86.2	±5	91.0	±5
	RDW-CV	%	16.0	±3.0	14.4	±3.0	15.0	±3.0
	HCT	%	17.6	±2.5	37.3	±3.5	45.8	±4.5
	PLT	×10 ⁹ /L	46	±20	258	±40	585	±60
	MPV	fL	8.5	±3.0	8.6	±3.0	8.6	±3.0